REBSTOCK®

Precision made in Germany

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GB

With the purchase of this instrument you obtain a high-quality product, whose appropriate application and handling are described below.

In order to keep dangers for patients and applicants to a minimum, it is very important to read this instruction manual carefully.

The application, disinfection, cleaning and sterilization must only be done by well-trained personnel.

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INSPECTION	<p>Before each use, the instruments must be carefully checked for their functionality. Instruments must not be used if there are surface damages like scratches, cracks, notches, kerfs, etc. as well as bent parts. Then the products have to be repaired or are to be given to the hospital waste disposal.</p> <p>NEVER USE ANY DAMAGED PRODUCTS!</p>
HANDLING	<p>Instruments must not be overstressed by distorting or levering, as this may lead to the damage or the breach of instrument parts.</p>
PURPOSE/ SCOPE OF APPLICATION	<p>Instruments as standard instruments for surgical interventions in general surgery.</p> <p>At any time, the surgeon and/or the applicant must decide if the instrument is suited for the specific application.</p> <p>The instruments must only be used by well-trained surgeons and/or applicants who are familiar with the necessary surgical technique of these medical products.</p> <p>Furthermore the surgeon is responsible for an appropriate training and sufficient information for the operating theatre staff as well as for an adequate expertise with the handling of the instruments.</p>
INDICATIONS	<p>No particular.</p>
CONTRAINDICATIONS	<p>None known.</p>
RISKS	<ul style="list-style-type: none">• Injuries of nerves, vessels, tissue• Bleedings• Infections• Thrombosis• Pulmonary embolism
COMPLICATIONS	<p>Generally complications occur rarely.</p> <p>The frequency and severity of complications depend on the kind of examination.</p>
COMBINATION WITH OTHER PRODUCTS/ INSTRUMENTS	<p>The products of REBSTOCK® INSTRUMENTS GmbH may in no case be combined with products, components and instruments of other manufacturers. Combinations with products of other manufacturers may have a negative impact on the result of the intervention and are not permitted, as the used components can possibly be not coordinated with each other. During the application it is recommended that only instruments and accessories of REBSTOCK® INSTRUMENTS GmbH are used.</p>
WASTE DISPOSAL	<p>In the case that instruments can no longer be repaired or re-processed, they should be given to the hospital-grade waste disposal.</p>
MATERIALS	<p>The raw materials used are stainless steels according to DIN EN ISO 7153-1, titanium alloys in accordance with DIN ISO 5832-3 or aluminium alloys in accordance with DIN EN 573-3.</p>

INSTRUCTION FOR REPROCESSING

TECHNIQUES	<ul style="list-style-type: none">• Cleaning• Disinfection• Sterilization with superheated steam (DIN EN ISO 17665-1)
WARNING	<p> The instruments are delivered in non-sterile condition and must be cleaned and, where appropriate, be disinfected and sterilized prior to use. Please also see our validated cleaning - and re-processing procedure!</p> <p> The instruments must only be re-processed by persons disposing of the necessary expertise and training. They must be able to assess the occurring risks with the related impacts.</p> <p> IN THE CASE THAT THE APPLICANT SHOULD DEVIATE FROM OUR PROCEDURE, THE SELECTED PROCEDURE HAS TO BE VALIDATED BY THE APPLICANT!</p>
LIMITATION OF REPROCESSING	<p>Frequent re-processing has a minor impact on the instruments. The end of the product lifetime is usually determined by wear and damages as a result of their use. The products must then be given to the hospital-grade waste disposal. Never use any damaged products!</p>

REBSTOCK® VALIDATED CLEANING- AND REPROCESSING PROCEDURE

PREPARATION FOR THE CLEANING AND DISINFECTION AT THE PLACE OF THEIR USE

- Remove visible surface contamination with a wet lint-free cloth.
- Flush non-visible surfaces like joints, hidden gaps, groves with distilled water, e. g. with a single-use syringe. It is recommended to re-process the instruments soonest possible after their use.
- Directly after use, a manual disinfection can be carried out in order to reduce the risk of an infection for the applicants. Though the instruments are inserted into a disinfectant liquid. It must be considered that the instruments are completely immersed in the disinfectant liquid, without blistering. Consider the manufacturer's instructions of the disinfectant liquid.

PREPARATION FOR THE DECONTAMINATION

Should instruments be detachable, detach them prior to their re-processing.

MANUAL CLEANING

Wherever possible, the cleaning and disinfection of instruments in cleaning devices should be preferred over the manual cleaning as automated procedures can be standardized!

- Immediately after use the instruments have to be cleaned in a mild alkaline cleaning liquid with a soft lint-free cloth, sponge and adequate brushes in order to remove blood and other residues.
- Do not use metal brushes or abrasive cleaners, as there is an increased risk of corrosion due to the damaged surface.
- The cleaning liquid used must be suited for the cleaning of steel-, titanium and aluminium products.
- The concentration and contact time of the cleaning liquid must strictly be followed according to the manufacturer instructions.
- Exchange the cleaning liquid at least once a day or immediately after a visible contamination.
- The cleaning liquid should be a non-foaming liquid, which is suited for manual cleaning.
- Disassemble and/or open the instruments as far as possible.
- Instruments with lumen and hollow spaces must be flushed carefully.
- Special attention should be paid to the cleaning of slots, catches, closures etc.
- Instruments with hinges and/or joints should be cleaned in open and closed position.
- Brushing and all further cleaning processes have to be carried out below the fluid surface in order to avoid the splashing of contaminated liquid.
- After the cleaning, the instruments have to be flushed intensively with clear running water. Remaining adhering residues have to be removed manually.
- The instruments should be stored in instrument cases which are suited for the cleaning, such as, for example, sterilization trays.

ULTRASONIC CLEANING

- For the support during the manual cleaning, for the removal of persistent residues before and after the automatic cleaning, for the cleaning of sensitive instruments, e. g. micro instruments.
- Use an ultrasonic device which is suited for the medical application.
- **THE INSTRUCTIONS OF THE DEVICE MANUFACTURER MUST BE FOLLOWED.**
- The wire basket of the ultrasonic device must be large and deep enough so that a complete immersion of the instruments is guaranteed.
- The instruments must be stored on a cleaning-friendly wire basket/sterilization tray in open position.
- Before cleaning, the ultrasonic cleaning bath should be heated up to a temperature, which is recommended from the manufacturers of the cleaning or cleaning-/disinfectant liquid.
- Generally, temperatures between 40°C and 50°C increase the cleaning effect.
- After the ultrasonic cleaning, the instruments have to be flushed with clear running water first. Demineralised water is recommended in the last rinse cycle. Afterwards the instruments have to be dried.
- The used cleaning-/disinfectant liquids must be prepared in accordance with the manufacturer's instructions.
- An overdosage of the cleaning-/disinfectant liquid must be avoided.
- The used cleaning-/disinfectant liquid must be suited for the cleaning of steel-, titanium and aluminium products.
- The cleaning liquid in the ultrasonic cleaning bath must be changed at least every working day. If is a visible contamination occurs several times a day.
- Products must be covered completely with the cleaning liquid.
- Use only wire baskets which do not affect the cleaning efficiency.
- Do not overload the sterilization trays.
- Prevention of „sonic shadows“.

Continued on next page

Continuation:

ULTRASONIC CLEANING

- All lumen and hollow spaces must be filled with cleaning liquid. It must be free of air bubbles.
- All lumen must be blown out with air in order to remove remainders of liquid.
- After the cleaning, the instruments must be macroscopically clean. That means free of visible contaminations.

MANUAL DISINFECTION

- The used disinfectant liquid must be suited for the disinfection of steel-, titanium and aluminium products.
- Cleaned instruments have to be placed in a tray with disinfectant liquid. The manufacturer's instructions concerning concentration, temperature and contact times have to be followed strictly.
- The disinfectant liquid must be freshly prepared every day. The multiple use can lead to an increased concentration caused by evaporation. This leads to a corrosion risk as well as to a high contamination level. The result is that the efficiency of the disinfection is reduced.
- Instruments with hinges and/or joints must be inserted in the disinfectant liquid in open position.
- All lumen and hollow spaces must be filled with disinfectant liquid. It must be free of air bubbles.
- Remove instruments from the disinfectant liquid with fresh single-use gloves.
- After the disinfection, the instruments have to be flushed sufficiently with clear running water.
- In order to avoid water spotting, the use of demineralised water is recommended.
- Instruments have to be dried directly after the rinse cycle.
- The instruments have to be free of visible contaminations after each cleaning process.
- Spotty instruments must be sorted out and require a special treatment (the cleaning procedure has to be repeated and/or an ultrasonic cleaning has to be carried out).
- Damaged or defective instruments must be sorted out immediately.

AUTOMATIC CLEANING AND DISINFECTION

Wherever possible, the cleaning and disinfection of instruments in cleaning devices should be preferred over the manual cleaning as automated procedures can be standardized!

- The specific characteristics of the device model must be considered.
- The used cleaning-/disinfectant liquid must be suited for the cleaning/disinfection of steel-, titanium and aluminium products.
- The cleaning liquid should be a non-foaming liquid.
- Instruments with problematic contaminations must be pre-treated manually or in the ultrasonic bath.
- Detach the instruments as far as possible prior to their reprocessing.
- For instruments with long lumen such as probes, an automatic cleaning is only suitable, if the hollow spaces can be flushed with the disinfection liquid and a secure rinsing can be guaranteed.
- The instrument cases like instrument trays should only be loaded to the extent that no negative impact on the cleaning process can occur.
- Load the instrument trays carefully so that all instruments can be flushed thoroughly (no rinse residues). The weight of the tray should not exceed 8 to 10 kg.
- After the removal from the cleaning device, the critical instruments with lumen, like cannulas, etc. must be analysed on possible remaining residues. The cleaning process has to be repeated in case of remaining residues and/or a manual cleaning has to be carried out (Please see the column „manual cleaning“).
- The machine-based thermal disinfection must be effected under consideration of the national requirements with regard to the A0-value (EN15883).
- After the end of the programme, it has to be checked if all programme levels were performed and if all control parameters are fulfilled.
- Open the cleaning-/disinfection device and remove the instruments with disinfected hands or with fresh single-use gloves.
- Dry lumen and hollow spaces with compressed air. Where necessary, dry the instruments with a lint-free cloth.

DRYING AND FUNCTIONAL CHECK

- Dry surfaces with a lint-free cloth as well as with compressed air.
- Dry all lumen and hollow spaces completely with compressed air.
- Check all instruments for their proper function.

**MAINTENANCE
AND CARE OF THE
INSTRUMENTS**

- The surface of the instruments must be dry and microscopically clean.
- Instrument oil (physiologically safe paraffin oil in accordance with DAB10 and/or PH.EUR or USP) must be applied on the surface, particularly on the moving parts/joints of the instruments.
- Sort out damaged instruments and send them in for repair.
- It is absolutely necessary that the instruments are sterilised before sending them in to the repair service.

PACKING

- Before sterilization, the instruments should be packed in a suited tray or a suited sterilization package (EN 868). The sterilization package is subject to the sterilization method, transport and storage. The package has a considerable influence on the sterilization result. The package has to be selected in that way that the instruments fit in the package. Sterilization packages have to be in accordance with the applicable standards.
- Use a sterilization indicator for each package and write down the sterilization- and expiry date on the package.

STERILIZATION

Recommended Sterilization Method	Pre-vacuum Stages	Recommended Temperature	Recommended Pressure	Holding Time	Drying Time
Steam pressure sterilization with a fractionated vacuum (DIN EN ISO 17665-1)	3 pre-vacuum stages with at least 60 millibar pressure	134 °C	3 bar	5 minutes	15 minutes

- During the sterilization, the instructions of the device manufacturer must be followed.
- After the sterilization, the sterile packaging and the sterilization indicators have to be checked on damages.
- The steam for sterilization purposes must correspond to DIN EN 285:2009. The steam which is used for the sterilization must be free of contaminations and must neither have a negative impact on the sterilization process nor cause any damage to the sterilizer or to the sterilized goods.

STORAGE

- Sterilized instruments have to be stored in a closed cabinet, protected from dust, humidity and temperature variations.

ADDITIONAL INFORMATION

Further indications for the re-processing of medical products:

www.rki.de

www.a-k-i.org

- ISO 17664, Sterilization of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices
- ISO 17665-1, Sterilization of health care products – moist heat, Part 1
- ISO 15883, Washer/Disinfectors: General Requirements, Terms and Definitions and Tests

**CONTACT WITH
THE MANUFACTURER**

Please see the Manufacturer- and Service Address

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The applicant is responsible that by means of well-trained personnel as well as by the application of suited devices and materials, the necessary re-processing result is achieved.

Guarantee

All products are manufactured from high-quality materials. Before dispatch every product passes our quality inspection. In case that there should occur a defect, please contact our service. The applicant has to decide about it by himself.

WE CANNOT ACCEPT ANY COMPENSATION CLAIMS IN CASE OF A NON-COMPLIANCE WITH THESE GUIDELINES AND/OR AN EVIDENT VIOLATION!

INAPPROPRIATE HANDLING AND USE OF OUR INSTRUMENTS CAN RESULT IN THE DAMAGE AS WELL AS THE FAILURE OF THE INSTRUMENT; THIS CAN CAUSE SERIOUS INJURIES OR EVEN THE DEATH OF THE PATIENT.

INAPPROPRIATE HANDLING OF THE SUPPLIED PRODUCTS WILL LEAD TO THE EXPIRY OF ANY GUARANTEE CLAIMS.

Caution

We do not assume any liability for the correct selection of the instrument which is suited for the surgical intervention.

In case that the instruments are used on patients with a Creutzfeldt-Jakob disease or a HIV infection, we disclaim any responsibility for the reprocessing.

REBSTOCK®
Precision made in Germany

REBSTOCK® INSTRUMENTS GmbH

In Weiheräcker 7
78589 Dürbheim-Tuttlingen
GERMANY

Fon +49 74 24 67 34

Fax +49 74 24 43 34

info@rebstock.de

/// www.rebstock.de